TABLE OF CONTENTS

English...	1	German...	9
French...	3	Product Codes...	11
Spanish...	5	Glossary of Symbols...	11
Italian...	7		

INTENDED USE

NPC-67[®] Neutralizing Buffer is used to neutralize the base pH reagents used in the N-acetyl-L-cysteine (NALC) digestion and decontamination procedure for the increased recovery of *Mycobacterium* spp. from sputum and other clinical specimens.

SUMMARY

The decontamination and digestion procedure, utilizing the compound N-acetyl-L-cysteine (NALC) combined with sodium hydroxide and sodium citrate (trisodium citrate), results in increased yields of tubercle bacilli. The NALC procedure utilizes N-acetyl-L-cysteine as a mucolytic compound by disrupting chemical bonds in mucus. The sodium hydroxide acts as a bacterial decontaminate and the sodium citrate (trisodium citrate) stabilizes the NALC by chelating (binding) any heavy metal ions present in the specimen. Since the sodium hydroxide has a pH of approximately 13.00, it will kill bacteria (including mycobacteria after 15–20 minutes of exposure). As such, timing of the decontamination is critical to limit the amount of *Mycobacterium* spp. killed by the basic pH. A pH indicator is incorporated into the decontamination reagent to monitor the pH throughout the decontamination and buffering procedure, allowing the laboratory technologist to visually see when neutralization has been achieved. Bringing the pH to a neutral range can stop the decontamination process. The NPC-67 is used to neutralize the NaOH following the appropriate digestion and decontamination time, resulting in a pH of below 8.10. Adding conventional M/15 Phosphate Buffer or phosphate buffered saline will result in a pH range of 9.40 to 12.20, requiring a titration to a neutral pH with 1N HCL, or continued decontamination of *Mycobacterium* spp. will occur. Studies have documented that pH values above 8.10 are toxic to *Mycobacterium* spp., including *Mycobacterium tuberculosis*. Following the decanting step, PRB[™] Pellet Resuspension Buffer is added to achieve a tight neutral pH value (6.80–7.10) in the specimen sediment, optimizing mycobacteria recovery.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

All clinical specimens submitted for the diagnosis of tuberculosis and other *Mycobacterium* spp. must be treated with appropriate care so as not to contaminate other specimens or laboratory personnel. Use all approved and regulated equipment for processing and detection procedures.

STABILITY AND STORAGE

NPC-67 is stable to the stated expiration date when stored at the required temperature. Prior to opening, NPC-67 can be stored at room temperature (15–30°C). After opening, store between 2–8°C. Allow the product to come to room temperature prior to use.

USER QUALITY CONTROL

Any product showing cloudiness, turbidity, precipitation or discoloration should be discarded. Quality controlled microorganisms should be utilized to verify procedures, media and reagents as appropriate for your laboratory's applicable regulatory agency or local procedural guidelines.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information.

PROCEDURE

Materials Provided: NPC-67 Neutralizing Buffer.

Materials Not Provided: Centrifuge, vortex mixer, sterile pipettes, microscope slides, TB media, centrifuge tubes, NAC-PAC[®] or NAC-PAC RED, PRB Pellet Resuspension Buffer, and CELL-BOND[®] Slides.

SPECIMEN PROCESSING

- Line up specimens (in centrifuge tubes) in a biosafety hood.
- Loosen specimen container caps. Work in sets equivalent to a centrifuge load.
- Open bottle labeled "NAC-PAC RED". Add the NALC powder to the NAC-PAC RED bottle. Shake well to dissolve the NALC powder. **NOTE:** Some residual NALC powder may remain in the vial. It is **not** necessary to liquefy the portion remaining in the vial. **THIS SOLUTION WILL BE GOOD FOR 72 HOURS AFTER MIXING.** Discard the mixed solution after 72 hours.
- To the sterile 50 ml centrifuge tube containing the specimen to be digested, add the NAC-PAC RED/NALC solution as follows:
 - For specimens **1–5 ml**, add a volume of NAC-PAC RED/NALC equal to that of the specimen volume.
 - For specimens **6–7 ml**, add 5 ml of NAC-PAC RED/NALC.
 - For specimens **8–10 ml**, add an equal volume of NAC-PAC RED/NALC and then split the specimen after step 6 equally into two centrifuge tubes, proceed with steps 7–9 and then combine the sediments from both tubes into one centrifuge tube and proceed with step 10.

Following this protocol will help achieve effective decontamination while also allowing for proper neutralization. If you routinely encounter specimens greater than 10 ml in volume, please contact Alpha-Tec Technical Services for special instructions.

- Tighten the caps on the centrifuge tubes. Mix each specimen on a vortex until liquefied (30 seconds per specimen).
- Allow each specimen to stand for 15–20 minutes. Vortex every 5 minutes during this step.
- Fill each tube with NPC-67 until effective neutralization is indicated by a color change from red/pink to colorless. Once a colorless point has been reached, do not continue to add NPC-67 to the sample. Tighten cap and swirl by hand to mix. **NOTE:** NPC-67 will achieve a colorless solution [basic pH neutralization] when added to a decontamination reagent [NAC-PAC RED] with a NaOH concentration of 3% or lower.
- Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended but not required to use a refrigerated centrifuge. Each laboratory must check the centrifuge head radius and use an appropriate nomogram for proper speed selection [rpm] to achieve the desired relative centrifugal field of 3000 xg.
- Working in a biosafety hood, pour off all supernatant into a splash-proof container holding an appropriate disinfectant. Use an appropriate disinfectant to disinfect any contamination on the lip of the specimen tube. Do not allow the disinfectant to run down inside the specimen tube.
- Resuspend the pellet with 0.5 ml–1.0 ml of PRB. Do not resuspend the pellet with NPC-67, water or saline. **NOTE:** To maximize time to detection for rapid growth automated detection systems, resuspend the pellet with 1.0 ml of PRB. Depending on the needs of your laboratory, the pellet may be resuspended with 0.5 ml of PRB to create a more concentrated sample for increased acid-fast smear staining sensitivity. Once the smears have been made, add an additional 1.0 ml of PRB to inoculate rapid broth detection systems and other media.
- Mix the sediment and buffer well, and inoculate the liquid broth for your automated detection equipment per the manufacturer's instructions.
- Place two drops of the sediment onto the surface of each of the TB media used. **NOTE:** A contamination control plate (BAP or TSA) can be inoculated at this point and incubated at 35–37°C for 48 hours.
- Make smears for acid-fast staining. Use adhesive CELL-BOND Slides or sterile albumin adhesive solutions to attach the specimen to the slide. Dry the smears and proceed with acid-fast staining per the manufacturer's directions. **NOTE:** An acid-fast stain control slide

should be stained in conjunction with the patient smears to verify the staining technique and components. Contact Alpha-Tec Technical Services for a complete list of acid-fast stains and control slides.

14. To the unused portion of the specimen, add the balance of the PRB and refrigerate at 2–8°C to save for future diagnostic procedures or reprocessing if necessary.

PROCEDURE NOTES

1. NPC-67 has been validated for use with multiple molecular diagnostic methods and systems. For more information regarding compatibility with specific methods or systems, contact Alpha-Tec Technical Services.
2. Small volume specimens with correspondingly low post-neutralization volumes can make centrifuge balancing difficult. If your laboratory frequently encounters small volume specimens, it is acceptable to add sterile saline to the sample to reach a combined volume of 5 ml prior to the addition of NAC-PAC RED/NALC solution. In this case, the sample should be decontaminated with 5 ml of NAC-PAC RED/NALC solution. This will increase the final post-neutralization specimen volume making centrifuge balancing easier.
3. Specimens contaminated with *Pseudomonas* spp. will need additional treatment with the OxA® Oxalic Acid Reagent Kit (#0004805). Refer to the Oxalic Acid Directions For Use for complete instructions, or call Alpha-Tec Technical Services for information on the pH effects of the oxalic acid procedure and the appropriate buffering requirements.
4. Following the decontamination of the specimen with NAC-PAC RED, bloody specimens may remain pink after the addition of the NPC-67 due to the residual hemoglobin in the specimen. If the color change cannot be visualized due to hemoglobin, add the NPC-67 up to the 50 ml mark to ensure complete neutralization. For additional information, contact Alpha-Tec Technical Services.

EXPECTED RESULTS

If *Mycobacterium* spp. are present in the clinical specimen and processed according to the procedures listed within this document, the recovery of cultivable, viable and clinically significant *Mycobacterium* spp. can be expected.

LIMITATIONS OF PROCEDURES

Timing of the decontamination step, proper buffering, speed and timing of the centrifugation step, proper decanting and addition of the PRB to the pellet are vital to the recovery of *Mycobacterium* spp. Failure to follow the listed procedures may result in decreased numbers or total loss of *Mycobacterium* spp. resulting in an inaccurate culture report.

SPECIFIC PERFORMANCE CHARACTERISTICS

NPC-67 was tested on clinical samples and recovered all culture appropriate *Mycobacterium* spp. when the designated procedures were followed.

BIBLIOGRAPHY

1. Babakhani, F., Warren N., Henderson D., and Dalton H. Effect of Transportation and Acid Neutralization on Recovery of Mycobacteria from Processed Specimens. *Am. J. Clin. Pathol.* 1995. 104(1):65–68.
2. Chapman J.S., Bernard, J.S. The Tolerance of Unclassified Mycobacteria. Limits of pH Tolerance. *Am. Rev Respir Dis.* 1962. 86:582–583.
3. Kent, P., Kubica, G.P., Public Health Mycobacteriology: A Guide for the Level III Laboratory. Centers for the Disease Control and Prevention, Atlanta. 1985.
4. Kubica, G.P., et al. Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria. *Am Rev Respir Dis.* 1963. 87:775–779.
5. Kubica, G.P., et al. Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of Mycobacteria. *Am Rev Respir Dis.* 1964. 89:284–286.

6. Lennette, E.H., et al., Manual of Clinical Microbiology. American Society for Microbiology Third Edition. 1980.
7. Vestal, A.L., Procedures for the Isolation and Identification of Mycobacteria. C.D.C., Atlanta, GA. 1975.
8. Yegian, D., Budd V. Toxic Effect of Sodium Hydroxide on Tubercle Bacilli. *Am. J. Clinical Pathology.* 1952. 22:456–460.

CONTACT

For Technical Assistance email Technical@AlphaTecSystems.com for Customer Service, email Sales@AlphaTecSystems.com or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8 am and 4 pm Monday through Friday, Pacific Time.

WARRANTY

CalibreScientific US, Inc. warrants this product to perform as described in the labeling and literature supplied. CalibreScientific US, Inc. disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall CalibreScientific US, Inc. be liable for any consequential damages arising out of aforesaid express warranty.

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BIBLIOGRAPHIE

1. Babakhani, F., Warren N., Henderson D., and Dalton H. Effect of Transportation and Acid Neutralization on Recovery of Mycobacteria from Processed Specimens. *Am. J. Clin. Pathol.* 1995. 104(1):65–68.
2. Chapman J.S., Bernard, J.S. The Tolerance of Unclassified Mycobacteria. Limits of pH Tolerance. *Am. Rev Respir Dis.* 1962. 86:582–583.
3. Kent, P., Kubica, G.P., Public Health Mycobacteriology: A Guide for the Level III Laboratory. Centers for the Disease Control and Prevention, Atlanta. 1985.
4. Kubica, G.P., et al. Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria. *Am Rev Respir Dis.* 1963. 87:775–779.
5. Kubica, G.P., et al. Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of Mycobacteria. *Am Rev Respir Dis.* 1964. 89:284–286.
6. Lennette, E.H., et al., Manual of Clinical Microbiology. American Society for Microbiology Third Edition. 1980.
7. Vestal, A.L., Procedures for the Isolation and Identification of Mycobacteria. C.D.C., Atlanta, GA. 1975.
8. Yegian, D., Budd V. Toxic Effect of Sodium Hydroxide on Tubercle Bacilli. *Am. J. Clinical Pathology.* 1952. 22:456–460.

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- Babakhani, F., Warren N., Henderson D., and Dalton H. Effect of Transportation and Acid Neutralization on Recovery of Mycobacteria from Processed Specimens. *Am. J. Clin. Pathol.* 1995. 104(1):65–68.
- Chapman J.S., Bernard, J.S. The Tolerance of Unclassified Mycobacteria. Limits of pH Tolerance. *Am. Rev Respir Dis.* 1962. 86:582–583.
- Kent, P., Kubica, G.P., Public Health Mycobacteriology: A Guide for the Level III Laboratory. Centers for the Disease Control and Prevention, Atlanta. 1985.
- Kubica, G.P., et al. Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria. *Am Rev Respir Dis.* 1963. 87:775–779.
- Kubica, G.P., et al. Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of Mycobacteria. *Am Rev Respir Dis.* 1964. 89:284–286.
- Lennette, E.H., et al., Manual of Clinical Microbiology. American Society for Microbiology Third Edition. 1980.
- Vestal, A.L., Procedures for the Isolation and Identification of Mycobacteria. C.D.C., Atlanta, GA. 1975.
- Yegian, D., Budd V. Toxic Effect of Sodium Hydroxide on Tubercle Bacilli. *Am. J. Clinical Pathology.* 1952. 22:456–460.

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BIBLIOGRAFIA

1. Babakhani, F., Warren N., Henderson D., and Dalton H. Effect of Transportation and Acid Neutralization on Recovery of Mycobacteria from Processed Specimens. *Am. J. Clin. Pathol.* 1995. 104(1):65–68.
2. Chapman J.S., Bernard, J.S. The Tolerance of Unclassified Mycobacteria. Limits of pH Tolerance. *Am. Rev Respir Dis.* 1962. 86:582–583.
3. Kent, P., Kubica, G.P., *Public Health Mycobacteriology: A Guide for the Level III Laboratory.* Centers for the Disease Control and Prevention, Atlanta. 1985.
4. Kubica, G.P., et al. Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria. *Am Rev Respir Dis.* 1963. 87:775–779.
5. Kubica, G.P., et al. Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of Mycobacteria. *Am Rev Respir Dis.* 1964. 89:284–286.
6. Lennette, E.H., et al., *Manual of Clinical Microbiology.* American Society for Microbiology Third Edition. 1980.
7. Vestal, A.L., *Procedures for the Isolation and Identification of Mycobacteria.* C.D.C., Atlanta, GA. 1975.
8. Yegian, D., Budd V. Toxic Effect of Sodium Hydroxide on Tubercle Bacilli. *Am. J. Clinical Pathology.* 1952. 22:456–460.

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LITERATUR

1. Babakhani, F., Warren N., Henderson D., and Dalton H. Effect of Transportation and Acid Neutralization on Recovery of Mycobacteria from Processed Specimens. Am. J. Clin. Pathol. 1995. 104(1):65–68.
2. Chapman J.S., Bernard, J.S. The Tolerance of Unclassified Mycobacteria. Limits of pH Tolerance. Am. Rev Respir Dis. 1962. 86:582–583.
3. Kent, P., Kubica, G.P., Public Health Mycobacteriology: A Guide for the Level III Laboratory. Centers for the Disease Control and Prevention, Atlanta. 1985.
4. Kubica, G.P., et al. Sputum Digestion and Decontamination with N-acetyl-L-cysteine-Sodium Hydroxide for Culture of Mycobacteria. Am Rev Respir Dis. 1963. 87:775–779.
5. Kubica, G.P., et al. Comments on the Use of the New Mucolytic Agent, N-acetyl-L-cysteine, as a Sputum Digestant for the Isolation of Mycobacteria. Am Rev Respir Dis. 1964. 89:284–286.
6. Lennette, E.H., et al., Manual of Clinical Microbiology. American Society for Microbiology Third Edition. 1980.
7. Vestal, A.L., Procedures for the Isolation and Identification of Mycobacteria. C.D.C., Atlanta, GA. 1975.
8. Yegian, D., Budd V. Toxic Effect of Sodium Hydroxide on Tubercle Bacilli. Am. J. Clinical Pathology. 1952. 22:456–460.

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Manufactured by CalibreScientific US, Inc.
1311 SE Cardinal Court, Suite 170
Vancouver, WA 98683 USA

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MDSS GmbH