TABLE OF CONTENTS

English...	1	German...	9
French...	3	Product Codes...	11
Spanish...	5	Glossary of Symbols...	11
Italian...	7		

INTENDED USE

Alpha-Tec Systems, Inc. OxA® Oxalic Acid Reagent Kit is used for decontamination of respiratory specimens containing *Pseudomonas* spp..

SUMMARY

Clinical specimens submitted to the laboratory for the isolation of acid-fast mycobacteria are often contaminated with commensal microbial flora. While most of the contaminating bacteria can be eliminated with a sodium hydroxide and N-acetyl-L-cysteine (NALC) procedure, *Pseudomonas* spp. can survive this decontamination procedure method, making accurate identification of mycobacteria extremely difficult. The use of oxalic acid as an effective decontamination method was first described by Corper and Uyei in 1930 — and later demonstrated by Whittier et al., in 1993, to decontaminate for *Pseudomonas aeruginosa* in specimens of cystic fibrosis patients.

FOR IN VITRO DIAGNOSTIC USE ONLY

PRECAUTIONS

- Oxalic Acid is POISONOUS. Avoid contact with skin and eyes. Should contact occur, flush immediately with water. Contact a physician if irritation occurs.
- NAC-PAC® RED contains a caustic chemical (sodium hydroxide). Use appropriate care in handling this reagent.
- Observe all laboratory standards for the handling of patient specimens. All clinical specimens submitted for the diagnosis of tuberculosis and other *Mycobacterium* spp. must be treated with the appropriate care to avoid contaminating other specimens or laboratory personnel. Use approved and regulated equipment for processing and detection procedures.

STABILITY AND STORAGE

Reagents in OxA are stable to the stated expiration date when stored at 15–30°C.

USER QUALITY CONTROL

Each bottle of reagent in OxA should contain the quantity listed on the label and should be clear and free of precipitates.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens for the detection of *Mycobacterium* spp. should be collected according to prescribed standards and delivered to the laboratory in a safe and timely manner. Refer to local procedural guidelines for this information.

PROCEDURE

Materials Provided: 5% Oxalic Acid, NAC-PAC RED, NPC-67® Neutralizing Buffer, PRB™ Pellet Resuspension Buffer.

Materials Not Provided: Centrifuge, centrifuge tubes (50 ml), microscope slides, sterile pipettes, vortex mixer, AFB Stains, TB Media.

SPECIMEN PROCESSING

ATTENTION: Handle all specimens with extreme caution. Possible Pathogens. Wear gloves at all times. Work in biosafety cabinet.

- Perform the NALC/NaOH digestion and decontamination procedure on the respiratory specimen. (Call your Alpha-Tec Account Executive for a complete list of our digestion/ decontamination and buffering products.)
- If patient has history of previous *Pseudomonas* spp. infection, or if the specimen is identified as having *Pseudomonas* spp.

contamination through a positive contamination control plate or rapid broth test, proceed with the OxA decontamination procedure.

OxA OXALIC ACID DECONTAMINATION PROCEDURE

- Use up to 3 ml from the original NALC/NaOH processed specimen and add it to a sterile 50 ml centrifuge tube.
- Add the entire contents of Oxalic Acid and vortex for 30 to 60 seconds.
- Allow to stand for 30 minutes, vortexing every 10 minutes.
- Add entire contents of NAC-PAC RED.
- After adding the NAC-PAC RED, mix gently.
- Add the entire contents of NPC-67. Tighten cap, and swirl to mix.
- Centrifuge the specimen tubes at 3000 xg for 15 minutes. It is recommended to use a refrigerated centrifuge. (Each laboratory must check the centrifuge head radius and use an appropriate nomogram for proper speed selection (rpm) to achieve the desired relative centrifugal field of 3000 xg.)
- Working in a biosafety hood, pour off all supernatant into a splash-proof container holding an appropriate disinfectant. Disinfect any contamination on the lip of the specimen tube. Do not allow the disinfectant to run down inside the specimen tube.
- Resuspend the pellet with 0.5–1.0 ml of PRB. Do not resuspend the pellet with any other reagents, water or saline.

Option 1: To maximize time to detection for rapid growth automated detection systems, resuspend the pellet with 1.0 ml of PRB.

Option 2: Depending on the needs of your laboratory, the pellet may be resuspended with 0.5 ml of PRB to create a more concentrated sample for acid-fast smear staining. Once the smears have been made, add an additional 1.0 ml of PRB for rapid broth detection methods and media inoculation.
- Mix the sediment and buffer well and inoculate the liquid broth for your automated detection equipment per the manufacturer's instructions.
- Place two drops of the sediment onto the surface of each of the TB media used. A contamination control plate (BAP or TSA) can be inoculated at this point and incubated at 35–37°C for 48 hours.
- Make smears for AFB staining. Use adhesive CELL-BOND® Slides or your appropriate sterile solutions to attach the specimen to the slide. Dry the smears and proceed with AFB staining per the manufacturer's directions. (Call your Alpha-Tec Systems Account Executive for a complete list of AFB Stains). **NOTE:** An AFB Quality Control Slide (Alpha-Tec #0003240) should be stained in conjunction with the patient smears to verify the staining technique and components.
- Add the balance of PRB to the unused portion of the specimen and refrigerate at 2–8°C for further diagnostic procedures or re-decontamination if the cultures or rapid detection methods indicate that a contaminate is present.

PROCEDURE NOTES

OxA Kit has been validated for use with multiple molecular diagnostic methods and systems. For more information regarding compatibility with specific methods or systems, contact Alpha-Tec Technical Services.

EXPECTED RESULTS

If Mycobacteria are present in the clinical specimen and processed according to the procedures listed within this document, the recovery of cultivable, viable and clinically significant Mycobacteria can be expected.

LIMITATIONS OF PROCEDURES

Timing of the decontamination step, proper buffering, speed and timing of the centrifugation step, proper decanting and addition of the PRB to the pellet are vital to the recovery of Mycobacteria. Failure to follow the listed procedures may result in decreased numbers of Mycobacteria or total loss of Mycobacteria resulting in an inaccurate culture report.

BIBLIOGRAPHY

- Corper, H.J., and N. Uyei. 1930. Oxalic Acid as a Reagent for Isolating Tubercule Bacilli and a Study of the Growth of Acid-Fast

- Nonpathogens on Different Mediums with Their Reactions to Chemical Reagents. J. Lab. Clin. Med. 15:348–369.
2. Lennette, E.H. et al., 1980. Manual of Clinical Microbiology. American Society of Microbiology, Washington, D.C., pp. 150–179.
 3. Vestal, A., 1975. Procedures for the Isolation and Identification of Mycobacteria. U.S. Dept. of Health, CDC Publication No. 79–8230, Center for Disease Control, Atlanta, GA. pp. 21–31.
 4. Whittier, S., et al., 1993. Improved Recovery of Mycobacteria from Respiratory Secretions of Patients with Cystic Fibrosis. J. Clin. Microbiol. 31:861–864.

CONTACT

Alpha-Tec Systems, Inc. offers a complete line of QC Slides, stains, reagents, and digestion systems for AFB specimen processing. For Technical Assistance email Technical@AlphaTecSystems.com and for Customer Service, email Sales@AlphaTecSystems.com or call [+1] 800.221.6058 or [+1] 360.260.2779 between 8 am and 4 pm Monday through Friday, Pacific Time.

WARRANTY

Alpha-Tec Systems, Inc. warrants this product to perform as described in the labeling and literature supplied. Any change or modification in the procedure may cause improper results. Unless this product is used in accordance with the labeling and literature, any and all warranties (expressed, implied or statutory) regarding the product including the implied warranty or merchantability and fitness for particular purpose are specifically denied. In no event shall Alpha-Tec Systems, Inc. be liable for any direct, indirect or consequential damages resulting from use or performance of this product.

TRADEMARKS

CELL-BOND®, NAC-PAC®, NPC-67®, OxA®, and PRB™ are trademarks of Alpha-Tec Systems, Inc., 1311 SE Cardinal Court, Suite 170, Vancouver, WA 98683 USA.

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BIBLIOGRAPHIE

1. Corper, H.J., and N. Uyei. 1930. Oxalic Acid as a Reagent for Isolating Tubercule Bacilli and a Study of the Growth of Acid-Fast Nonpathogens on Different Mediums with Their Reactions to Chemical Reagents. *J. Lab. Clin. Med.* 15:348–369.
2. Lennette, E.H. et al., 1980. *Manual of Clinical Microbiology*. American Society of Microbiology, Washington, D.C., pp. 150–179.
3. Vestal, A., 1975. *Procedures for the Isolation and Identification of Mycobacteria*. U.S. Dept. of Health, CDC Publication No. 79–8230, Center for Disease Control, Atlanta, GA. pp. 21–31.
4. Whittier, S., et al., 1993. Improved Recovery of Mycobacteria from Respiratory Secretions of Patients with Cystic Fibrosis. *J. Clin. Microbiol.* 31:861–864.

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BIBLIOGRAFIA

1. Corper, H.J., and N. Uyei. 1930. Oxalic Acid as a Reagent for Isolating Tubercule Bacilli and a Study of the Growth of Acid-Fast Nonpathogens on Different Media with Their Reactions to Chemical Reagents. *J. Lab. Clin. Med.* 15:348–369.
2. Lennette, E.H. et al., 1980. *Manual of Clinical Microbiology*. American Society of Microbiology, Washington, D.C., pp. 150–179.
3. Vestal, A., 1975. *Procedures for the Isolation and Identification of Mycobacteria*. U.S. Dept. of Health, CDC Publication No. 79–8230, Center for Disease Control, Atlanta, GA. pp. 21–31.
4. Whittier, S., et al., 1993. Improved Recovery of Mycobacteria from Respiratory Secretions of Patients with Cystic Fibrosis. *J. Clin. Microbiol.* 31:861–864.

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Manufactured by Alpha-Tec Systems, Inc.
MDSS GmbH